

B337 Instant Wide LTC Low Bed

Bariatric Hi-Lo Bed with Width and Length Expansion


The B337 is one of the most versatile bariatric Hi-Lo Beds available. With built-in toolless length and width adjustment and 600 lb. capacity the B337 adapts to every need.



ADJUSTABLE WIDTH



ADJUSTABLE LENGTH



EMBEDDED FOOTBOARD CONTROL AND HAND PENDANT



ONE-STEP LOCK & UNLOCK

- Whisper quiet actuators move the deck from a low height of just 7", to a safe working height of 30".
- The sleep surface easily expandable from a standard 36" width to 39" and 42" widths.
- Length is also adjustable from a standard 80" to 84" and 88" for the tallest residents.
- The B337 weight capacity is 600 lbs.
- 4 Function bedside hand pendant with backlighting.
- The staff control panel on the footboard includes all pendant functions with the addition of Trendelenberg, reverse Trendelenberg, underbed lights, lock-out, auto-contour and CPR.
- The B337 easily rolls at any height and is easily locked or unlocked with color coded foot pedals.
- Standard features include: adj. mattress retainers, heel lift, wall guard, IV pole holders
- Assist rails and swing away half-rails available.
- The B337 is delivered standard with a 36" or 42" foam core mattress.
- Streamline Medequip can also provide additional bariatric and patient management products such as: Bariatric Air Mattress Systems, Heavy Duty Patient Lifts, Bariatric Wheelchairs and Fall Pads.



(346) 445-4017

www.streamlinemedequip.com