


Protekt Protopia™ EXP

Expandable Ultra-Low Long Term Care Bed



9" Low bed height



Tool free width expansion



Tool free length extension



10 Function hand control

- Lowest position of only 9"
- Heavy duty welded steel construction.
- 10 function hand pendant.
- Trendelenburg & reverse Trendelenburg function.
- Auto contour function
- Adjustable extending leg elevator
- Pair of quarter rails or assist bars and head & foot boards included.
- Integrated tool free width and length extension.
- Bumper bar protects walls from damage.
- Corner mattress retainers.
- Quick release head and foot boards.
- (8) 3" casters: 4 locking and 4 swivel.
- Weight Capacity: 750 lbs.
 - Warranty:
 - 5 Years on Electrical Components
 - 15 Years on Frame (1 Year on all other parts)



Quarter Rail



Pivot Assist



Q-Bar



Soft Touch Rotating



Mahogany



Amber Cherry



Oak

Streamline
MedEquip

www.streamlinemedequip.com

346-445-4017