


Padded Divided Leg Sling w/o Head Support



| ITEM # | DESCRIPTION |
|------------------------------------|---|
| 30100P-NHS Small | Solid padded divided leg w/o head support: small 600 lb. weight capacity |
| 30100-NHS Medium | Solid padded divided leg w/o head support: medium 600 lb. weight capacity |
| 30101-NHS Large | Solid padded divided leg w/o head support: large 600 lb. weight capacity |
| 30102-NHS X-Large | Solid padded divided leg w/o head support: X-large 600 lb. weight capacity |

*All slings have been independently weight capacity certified by Intertek Testing Services according to ISO10535 standards.

