

PROTEKT® 500 STAND

34500 | SIT-TO-STAND PATIENT LIFT

The Protekt® 500 lb. Sit-to-Stand patient lift has a higher weight capacity than most other competitors sit-to-stand lifts. With its 3" dual front casters it ensures safe and easy transfers. Our performance class actuator includes an emergency quick release.

HCPCS
E0637

- 4 sling hooks allows for multiple sling options, safer patient handling and a wider range of transfers.
- Fast speed lift motor reduces time for the transfer.
- Battery lasts twice as long as competing brands for less downtime (24 volts).
- Ergonomically designed foot pedal offers hands free opening of base.
- 4" rear casters with brake & 3" dual front casters.
- Single piece molded kneepad with handlebar.
- Height & depth adjustable kneepad position.
- Full parts breakdowns and modular design on lift and electronics.
- Advanced electronics with auxiliary up/down lift controls.

WARRANTY:
Electronics 2 years
Frame 3 years
Battery 1 year

SPECIFICATION FOR ITEM # 34500

Weight Capacity	500 lbs.
Unit Weight	119 lbs.
Base Ext. Width - Closed	26.4"
Mast Height	41.7"
Base Length	41.3"
Front Caster Height	4.7"
Knee Pad Depth Range	2.8"
Knee Pad Height Range	20.2"-23.7"
Footplate Width	13.8"
Footplate Length	14.5"
Base Opened Turning Diameter	54.2"
Base Closed Turning Diameter	50.0"

Battery	LP-TC12-BPAK
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN

