

Protekt® Aire 6000

Low Air Loss/Alternating Pressure Mattress System
w/ Deluxe Digital Pump, Cell-On-Cell Support Base



HCPCS
E0277

FEATURES AND BENEFITS:

PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating cycles 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Auto firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

80060	Protekt® Aire 6000 36"x80"x8"
80065	Protekt® Aire 6000 with Raised Rails

MATTRESS & COVER

- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Available in 84" length.
- Optional raised side rails.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

