

Protekt Protopia™

Ultra-Low Long Term Care Bed



7" Low bed height

Best in class in several function categories*

- *Lowest position of only 7" (7" to 31" range)
- Heavy duty welded steel construction.
- *8 function hand pendant.
- *Auto contour function.
- Pair of quarter rails or assist bars and head & foot boards are included.
- *Integrated tool free length extension. (80"- 84"- 88")
- Bumper bar protects walls from damage.
- Corner mattress retainers.
- Quick release head and foot boards.
- (8) 3" casters: 4 locking and 4 swivel.
- *Weight Capacity: 500 lbs.
- *Warranty: 5 Years on Electrical Components, 15 Years on Frame
- *Bed end and choice of rails included in price



Tool free length extension



8 Function hand control



Quarter Rail



Pivot Assist



Q-Bar



Soft Touch Rotating



Mahogany



Amber Cherry



Oak

Streamline
MedEquip

www.streamlinemedequip.com

346-445-4017