

# Full Body Mesh Slings



ITEM #	DESCRIPTION
<b>30110</b> <b>Medium</b>	Mesh full body: medium 600 lb. weight capacity
<b>30111</b> <b>Large</b>	Mesh full body: large 600 lb. weight capacity
<b>30124</b> <b>X-Large</b>	Mesh full body: X-large 600 lb. weight capacity
<b>30118</b> <b>XX-Large</b>	Mesh full body: XX-large 600 lb. weight capacity

\*All slings have been independently weight capacity certified by Intertek Testing Services according to ISO10535 standards.

