


# Protekt Protopia™ EXP

## Expandable Ultra-Low Long Term Care Bed



7" Low bed height



Tool free width expansion



Tool free length extension



8 Function hand control

### Best in class in several function categories\*

- \*Lowest position of only 7" (7" to 31" range)
- Integrated tool free width expansion. (36", 42", 48")
- \*Integrated tool free length extension. (80"- 84"- 88")
- Heavy duty welded steel construction.
- \*8 function hand pendant.
- \*Auto contour function.
- Bumper bar protects walls from damage.
- Corner mattress retainers.
- Quick release head and foot boards.
- (8) 3" casters: 4 locking and 4 swivel.
- \*Weight Capacity: 550 lbs.
- \*Warranty: 5 Years on Electrical Components, 15 Years on Frame
- \*Bed ends and choice of rails included in price (Soft touch rails extra)



Quarter Rail



Pivot Assist



Q-Bar



Soft Touch Rotating



Mahogany



Amber Cherry



Oak

**Streamline**  
MedEquip

[www.streamlinemedequip.com](http://www.streamlinemedequip.com)

**346-445-4017**